



REORDER NUMBERS:

10165 – WITHOUT WHEELS

10166 – WITH WHEELS

BARIATRIC DUAL RELEASE FOLDING WALKER

500 LB. WEIGHT CAPACITY

- Sturdy aluminum structure provides reliable stability, while remaining lightweight
- Tool-free push button allows for height adjustments in 1-inch increments
- Contoured handgrips allow for added patient comfort
- Opens and closes with a push-button mechanism, easily operated by fingers or palms
- Two-button folding capabilities allow patients to fold each side independently
- 10166: Dual 5-inch front wheels provide increased ease of movement on a variety of surfaces
- Limited Lifetime Warranty
- HCPCS: E0149



10165



10166