



REORDER NUMBER: 10169



## BARIATRIC HD STEEL WALKER

700 LB. WEIGHT CAPACITY

- Dual-purpose walker can be customized as standard walker or rolling walker
- Comes equipped with 4 stabilizing legs and 2 legs with wheels
- Sturdy steel structure provides reliable stability
- Strong front cross braces and dual side braces provide added support
- Designed with an extra-wide frame for a more comfortable fit
- Tool-free push button allows for height adjustments in 1-inch increments
- Contoured handgrips allow for added patient comfort
- Opens and closes with a push button mechanism, easily operated by fingers or palms
- Two-button folding capabilities allow patients to fold each side independently
- Limited lifetime warranty

