



Folding Bedside Commode

The DME Care Solutions Folding Bedside Commode is a three-in-one unit that functions as a commode, raised toilet seat, or toilet safety frame. Designed with a snap-on seat, tool-free removable back, and drop-in commode bucket with splash guard, it folds flat for convenient storage and transport.

- Versatile, three-in-one commode opens easily and folds flat
- Constructed from powder-coated steel for optimal durability
- Can be configured as a commode, raised toilet seat, or toilet safety frame over existing toilets
- Includes durable plastic snap-on seat with lid, drop-in bucket with handle, splash guard, plastic armrests for added comfort, and removable back
- 10372C and 10373C: Feature collapsible commode bucket for compact storage and portability
- Tool-free setup and 1-inch height adjustments (from 16.22" to 22.24") accommodate a range of users
- Weight capacity 350 lb.
- Comes in retail-friendly packaging
- Limited Lifetime Warranty

Item	Description	Qty.	HCPC
10370	Standard, Bulk Packed	4 / box	E0163
10371	Elongated, Bulk Packed	4 / box	E0163
10372C	Standard, Collapsible Bucket	1 / box	E0163
10373C	Elongated, Collapsible Bucket	1 / box	E0163