



REORDER NUMBER: 10524



BARIATRIC GERI CHAIR RECLINER

3-POSITION • BLUERIDGE | 500 LB. WEIGHT CAPACITY

- Neutral-colored, sturdy steel frame provides optimal support
- Three chair positions for optimal flexibility
- Retractable lock bar secures the chair in the desired position
- Wide, padded armrests and built-in headrest for maximum support and comfort
- Blow-molded tray locks in 5 positions
- Side panels are removable for easy cleaning and maintenance
- 5" casters (2 locking and 2 non-locking)
- Moisture barrier on vinyl seat helps prevent fluid seepage
- Foam and vinyl meet California Technical Bulletin #117
- Warranty: Limited lifetime on frame; 1 year on upholstery and moving parts
- HCPCS: E1031

