



L-Mesitran® Tulle Honey Dressing

Impregnated with
L-Mesitran Soft Gel
Containing Medical-Grade Honey



ITEM 3071 4" × 4" (10 cm × 10 cm)

Dürma+ L-Mesitran Tulle is a non-adhering mesh dressing impregnated with a honey-based wound gel that helps manage exudate while maintaining a moist healing environment. The medical-grade honey formula provides antimicrobial support and helps manage odor.

- Non-adhering polyethylene mesh dressing provides gentle wound contact
- Honey-based gel helps to draw exudate away from wound bed to maintain balance
- Infused with 40% medical-grade honey for broad-spectrum antimicrobial support
- Gel composition supports autolytic debridement
- Odor-control properties assist with management of fungating wounds
- Blend of honey, hypo-allergenic lanolin, and vitamins C & E helps maintain moisture
- Designed for use with secondary dressings to reduce discomfort during changes
- Suitable for chronic and acute wounds, burns, ulcers, and donor sites
- 10 per box / 10 boxes per case



| 1.888 DYNAREX | DYNAREX.COM |

| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 |