



SINGLE SUCTION CATHETERS

Graduated with Control Valve



Dynarex Resp-O₂™ Single Suction Catheters are designed to safely and effectively remove secretions from the respiratory tract of tracheostomy patients. The soft, flexible Catheter with beveled tip is gentle on delicate mucosal tissue yet firm enough to allow for catheter rotation. The Catheter includes a vented cap port for enhanced suction control and two staggered eyelet openings to increase suctioning effectiveness and reduce patient trauma.

- Soft, pliable catheter with beveled tip provides optimal patient comfort
- Vented cap with control valve allows for greater procedural control
- Catheter tip is designed with 2 staggered eyelet openings for maximum suctioning efficiency
- Catheter features easy-to-read depth-number markings
- Sterile, single-patient use
- Not made with natural rubber latex or DEHP



6 Fr

Pediatric

COIL PACKED

ITEM 32000

STRAIGHT PACKED

ITEM 32000S

PACKAGE QUANTITY
50/Case

8 Fr

Pediatric

COIL PACKED

ITEM 32001

STRAIGHT PACKED

ITEM 32001S

PACKAGE QUANTITY
50/Case

10 Fr

Pediatric / Adult

COIL PACKED

ITEM 32002

STRAIGHT PACKED

ITEM 32002S

PACKAGE QUANTITY
50/Case

12 Fr

Adult

COIL PACKED

ITEM 32003

STRAIGHT PACKED

ITEM 32003S

PACKAGE QUANTITY
50/Case

14 Fr

Adult

COIL PACKED

ITEM 32004

STRAIGHT PACKED

ITEM 32004S

PACKAGE QUANTITY
50/Case

16 Fr

Adult

COIL PACKED

ITEM 32005

STRAIGHT PACKED

ITEM 32005S

PACKAGE QUANTITY
50/Case

18 Fr

Adult

COIL PACKED

ITEM 32006

STRAIGHT PACKED

ITEM 32006S

PACKAGE QUANTITY
50/Case