



PORTABLE OXYGEN CONCENTRATOR

ITEM 33955

OXYGEN CONCENTRATION
90% (-3% +6%)

PULSE DOSE RANGE
1-6 L/Min

DECIBEL LEVEL
43 dBA

BATTERY LIFE
Up to 4 hours on flow setting 1

WEIGHT WITHOUT BATTERY
4.01 lb.

WEIGHT WITH BATTERY
5.05 lb.

PACKAGE QUANTITY
1/Case

The Resp-O₂™ Portable Oxygen Concentrator is a premium respiratory device designed to effectively administer supplemental oxygen. It delivers 90% oxygen concentration at flow rates of 1-6 liters per minute to suit a range of needs. The FAA-compliant concentrator features a compact, lightweight design for convenient travel and storage.

- Portable concentrator designed to provide consistent supplemental oxygen at 90% concentration
- Adjustable flow settings of 1-6 liters per minute accommodate a range of therapeutic needs
- Small footprint allows use in vehicles, offices, and bedside spaces
- Lightweight construction allows for convenient transport, travel, and storage
- FAA-compliant design supports air travel and long-distance mobility
- Visual and audio alarms alert users to power disruptions, low oxygen purity, and lack of breath detection
- Simple, intuitive controls allow for easy operation
- Typical sound level of 43 dBA (maximum 66 dBA) helps reduce disruption during use
- AC/DC adapter and high-capacity battery provide flexible power options
- Carrying bag with shoulder strap supports hands-free portability



INCLUDED:

Carrying Case, Shoulder Strap, Power Cables, and Tubing

