



DEAD END CAPS

Luer Lock

ITEM 7057

Dynarex Dead End Caps offer a safe and effective solution for sealing dialysis ports after treatment. These sterile, non-vented caps help minimize the risk of contamination by securely closing the blue and red ports without an access port. Each cap is designed with wings to enhance grip and allow for precise handling.

- Sterile caps designed to close blue and red ports after dialysis
- Non-vented design without access port helps to reduce the risk of infection
- Crafted with wings for improved grip and control
- Not made with natural rubber latex or DEHP
- Packaged as 2 units per pouch, 1 red cap and 1 blue cap
- 100 pouches per box, 10 boxes per case

