

THE MONTCLAIR COLLECTION



LONG TERM CARE BED RAILS & BARS

- Assistive bars and rails provide stability to users getting in and out of bed
- Compatible with all Dynarex D-Series LTC beds
- **Rotating Assist Bar:** Durable steel construction that rotates 90° in either direction
- **Pivot Assist Bar:** Rotates down when not in use for added freedom of movement
- **Metal Swing Rail:** Easy-to-use push-button mechanism allows rail to swing up and down
- **Composite Swing Rail:** Pleasing design and comfort-touch technology provide a gentle feel
- Universal Bar options available for easy attachment to most brands' LTC beds
- Meets entrapment guidelines when properly installed
- 1 pair per case

12950	METAL SWING RAIL	1 PAIR/CS
12951	COMPOSITE SWING RAIL	1 PAIR/CS
12952	ROTATING ASSIST BAR	1 PAIR/CS
12953	PIVOT ASSIST BAR	1 PAIR/CS
12960	UNIVERSAL ROTATING ASSIST BAR	1 PAIR/CS
12961	UNIVERSAL PIVOT ASSIST BAR	1 PAIR/CS