



| 1.888 DYNAREX | DYNAREX.COM |

COMBO PACK

## ACTISPLINT WITH SENSI-WRAP ROLL

BREATHABLE, RADIOLUCENT

ITEM **3533**

- Conformable splints are designed to help stabilize fractures
- Aluminum core allows for optimal support and durability
- Breathable splints feature outer foam layers for improved comfort
- Radiolucent materials allow for unobstructed patient imaging
- Easy to cut with scissors for added customization
- Features water-resistant design for increased protection
- Fastens easily with most tapes or wraps
- Lightweight splints are ideal for first-aid, emergencies, sports, and more
- CAUTION: This product contains natural rubber latex which may cause allergic reactions
- DEHP Free
- 50 per case



### INCLUDES

- 36" Rolled ActiSplint (Red/Charcoal Gray)
- 2" x 5 yd. Self-Adherent Bandage Roll (Black Sensi-Wrap)



| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 |