



REORDER NUMBERS:


10431 – PLUS, 42"×80"

10432 – PLUS, 48"×80"

10437 – HD MULTIZONE, 54"×80"

BARIATRIC FOAM MATTRESSES

- Designed to help prevent and manage pressure wounds
- 10431-10432: Constructed from three layers of variable density foam
- 10431-10432: Foam top features precision diamond cuts and memory foam heel zone
- 10437: Constructed with variable density foam for optimal support
- Two-way stretch cover is designed to reduce shear friction and facilitate air flow
- Cover is fluid-resistant with anti-contamination zipper for optimal hygiene
- 42", 48", and 54" widths help to accommodate larger patients
- Limited Warranty: 3-year on Foam, 1-year on Cover
- HCPCS: E0184



WEIGHT CAPACITY

- **600 lb.** (10431, 10432)
- **1000 lb.** (10437)