

SafeTouch™

BLACK NITRILE EXAM GLOVES

Item 2521 – Small

Item 2522 – Medium

Item 2523 – Large

Item 2524 – X-Large

- **Chemo & Fentanyl Tested per ASTM D6978-05**
- High-quality nitrile gloves crafted for optimal protection and comfort
- Designed for tactile sensitivity and puncture-resistance
- Tear-resistant polymer lining provides added durability
- Textured fingertips provide improved grip for enhanced performance
- Thickness range: 3-4 mils at palms and 4-6 mils at fingertips
- Secure fit around wrist and hand-conforming flexibility help prevent fatigue
- Beaded cuffs allow for easy donning
- Ambidextrous and powder-free
- Ideal for tattoo studios, law enforcement, medical settings, EMS, beauty salons, and more
- Not made with natural rubber latex
- 100 per box, 10 boxes per case