




PEDIATRIC CLOSED SUCTION CATHETERS

Resp-O₂™ Closed Suction Catheters are tailored for neonatal and pediatric patients, offering safe and effective suctioning while maintaining continuous ventilation and minimizing infection risks.

- Designed for neonatal and pediatric patients to support delicate airway care
- Closed suction system maintains sterility and continuous ventilation, reducing infection risk
- Color-coded depth markers provide visual cues for safe and accurate catheter placement
- Smooth, flexible tubing allows for more comfortable insertion and greater ease-of-use
- Includes lavage/rinse ports and lockable thumb control valve for convenient cleaning and added protection
- Not made with DEHP



ACCESSORIES

ELBOW

Suction Depth Stickers

Change of Day Sticker

31062: 3.0, 3.5 and 4.0 Endo Tube Adaptors

Y ADAPTOR

Suction Depth Stickers

Change of Day Sticker

3.0, 3.5 and 4.0 Y Adaptors

8 Fr	10 Fr	8 Fr	10 Fr
TRACHEAL	TRACHEAL	ENDOTRACHEAL	ENDOTRACHEAL
Elbow • 20.3cm / 8"	Elbow • 20.3cm / 8"	Y Adaptor • 30.5cm / 12"	Elbow • 40.5cm / 16"
ITEM 31053	ITEM 31054	ITEM 31061	ITEM 31063
		Elbow • 35.5cm / 14"	
		ITEM 31062	
PACKAGE QUANTITY 20/Case	PACKAGE QUANTITY 20/Case	PACKAGE QUANTITY 20/Case	PACKAGE QUANTITY 20/Case