



DRAPE SHEETS

Dynarex Drape Sheets are designed to provide patients with personal privacy during medical exams. These drape sheets feature 2-ply paper construction for enhanced strength and are available in a variety of sizes and colors.

- Designed to provide patients with personal privacy during exams and procedures
- Constructed from 2-ply paper for added strength and reliability
- Soft and breathable tissue offers improved comfort
- Single-use, disposable sheets allow for optimal hygiene
- Available in a variety of sizes and attractive colors
- 8150: Features polyethylene layer for added protection against fluids
- Not made with natural rubber latex

8121	40" × 48"	White	100 / cs
8122	40" × 48"	Blue	100 / cs
8123	40" × 48"	Mauve	100 / cs
8131	40" × 60"	White	100 / cs
8132	40" × 60"	Blue	100 / cs
8133	40" × 60"	Mauve	100 / cs
8141	40" × 72"	White	50 / cs
8150	40" × 90"	Blue	50 / cs
8150BLK	40" × 90"	Black	50 / cs
8151	40" × 90"	White	50 / cs
8156	40" × 90"	Black	50 / cs

