



IV ADMIN AND EXTENSION SETS

- Available in a variety of lengths and configurations for optimal adaptability
- For gravity infusion only; not compatible with infusion pumps EXCEPT Item # 7068 which is compatible with Zyno Pump
- Sterile and individually packaged for optimal hygiene
- Non-pyrogenic; Not made with natural rubber latex or DEHP

ADMINISTRATION SETS

- Designed to help deliver fluids and medications safely and effectively
- Connectors can sustain water pressure of 200 kpa for 15 minutes and air pressure of 50 kpa for 15 seconds
- 50 / bx (7029: 40 / bx)

EXTENSION SETS

- Constructed from premium plastic materials for improved reliability
- Designed to add length to IV delivery systems
- Tubing can sustain static tensile force of 15 N for 15 seconds
- 100 / bx



IV ADMINISTRATION SETS

10 drops 15 drops 20 drops 60 drops

Item	Drops / mL	Length	Needleless Port	Pre-Pierced Port	Vented	Micron Filter	Flow Regulator	Accurate Regulator	Notes	Illustration
7029	10	98" (249 cm)	2		✓		✓			
7030	10	83" (211 cm)	1	1						
7031	15	83" (211 cm)	1	1						
7032	15	78" (198 cm)		1						
7033	15	106" (269 cm)		2						
7034	15	73" (185 cm)		1						
7037	15	92" (234 cm)		2						
7041	10	89" (226 cm)		2					Primary Set with Secondary / Piggyback Port and Check Valve	
7042	10	112" (284 cm)	2						Primary Set with Secondary / Piggyback Port and Check Valve	
7044	20	92" (234 cm)		1	✓					
7045	20	83" (211 cm)	1		✓	0.2		✓		
7048	20	100" (254 cm)		1		0.2				
7049	20	92" (234 cm)		1	✓	0.2	✓			
7051	60	83" (211 cm)	1	1						
7052	60	78" (198 cm)		2	✓					
7055	20	84" (213 cm)	1			15		✓		
7068	20	106" (269 cm)							Compatible with Zyno Pump	

IV EXTENSION SETS (NEEDLE-FREE LUER LOCK CONNECTOR)

Item	Length
7060	6" (15.24 cm)

Item	Length
7061	7" (17.78 cm)

Item	Length
7062	8" (20.32 cm)



| 1.888 DYNAREX | DYNAREX.COM |

| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 |