



OVERBED TABLES

ITEM 10450-10456

The Dynarex Overbed Table is a durable, wheeled table designed to provide reliable support and convenient access to meals and personal items. The table features a height-adjustable base, laminate wood tabletop, and locking wheels for controlled movement and stability.

- Provides reliable bedside support and access to meals and personal items
- Steel-welded H-base frame with powder-coated finish provides stable support
- Table height adjusts from 28" to 45" to fit most standard-height beds
- Spring-loaded locking lift mechanism allows easy adjustment and secure placement
- Locking wheels allow smooth movement and help keep the table in place
- Laminate wood tabletop provides a durable, easy-to-clean surface
- Supports up to 50 lb. distributed weight capacity
- Ideal for long-term care, rehabilitation, assisted living, and homecare settings
- HCPC: E0274





OVERBED TABLES

10450-WN	Walnut	1 / cs
10450-MH	Mahogany	1 / cs
10450-OA	Oak	1 / cs
10450-CH	Cherry	1 / cs

LOW OVERBED TABLE

10456	Walnut	1 / cs
-------	--------	--------

