

SafeTouch™

BLUE NITRILE EXAM GLOVES

Item 2511 – Small

Item 2512 – Medium

Item 2513 – Large

Item 2514 – X-Large



- **Chemo & Fentanyl Tested per ASTM D6978-05**
- High-quality nitrile gloves crafted for optimal protection and comfort
- Designed for improved tactile sensitivity and puncture-resistance
- Tear-resistant polymer lining provides added durability
- Textured fingertips provide improved grip for enhanced performance
- Thickness range of 3-4 mil at palm, 4-6 mil at fingertips
- Secure fit around wrist and hand-conforming flexibility help prevent fatigue
- Beaded cuffs allow for easy donning
- Ambidextrous and powder-free
- Not made with natural rubber latex
- 100 per box, 10 boxes per case