

*Safe***Touch™**

STERILE NITRILE SURGICAL GLOVES

Item 6537 – 8.0

Item 6538 – 8.5

Item 6539 – 9.0



- High-quality sterile nitrile gloves crafted for optimal protection and comfort
- Designed for improved tactile sensitivity and puncture-resistance
- Tear-resistant polymer lining provides added durability
- Thickness range of 5-8 mil at the palm and fingertips
- Secure fit around wrist and hand-conforming flexibility help prevent fatigue
- Beaded cuffs allow for easy donning
- Gloves are hand-specific and powder-free
- Not made with natural rubber latex
- 50 pairs per box, 4 boxes per case