



10261
HCPCS A4358



10259
HCPCS E0978

10262
HCPCS E0978



36622
HCPCS E1354



10559
HCPCS E0776



10260
HCPCS E0971



10269
HCPCS E0961



10270
HCPCS E0961

BARIATRIC WHEELCHAIR ACCESSORY COMPATIBILITY CHART

	STANDARD 10233-10234	HD 10290-10293	RECLINING 15015-24	RECLINING 15015-26	RECLINING 15015-28	TRANSPORT 10244
10260 ANTI-TIPPER	YES	YES	NO	NO	NO	YES
10261 URINARY DRAINAGE BAG HOLDER	YES	YES	YES	YES	YES	YES
10259 SEATBELT, 72"	YES	YES	YES	NO	NO	YES
10262 SEATBELT, 60"	YES	YES	YES	NO	NO	YES
10559 TELESCOPING IV POLE, 40"-75"	YES	NO	NO	NO	NO	YES
10269 BRAKE LEVER EXTENSION, 8"	YES	YES	YES	NO	NO	YES
10270 BRAKE LEVER EXTENSION, 11"	YES	YES	YES	NO	NO	YES
36622 OXYGEN TANK HOLDER	YES	NO	NO	NO	NO	YES

SEATBELTS

- Made from durable and comfortable fabric with adjustable lock-and-release straps

BRAKE LEVER EXTENSIONS

- Easily slips onto and adds length to existing standard wheel lock lever

IV POLE

- Durable steel construction
- Easy to install and assemble

OXYGEN TANK HOLDER

- Oxygen tank holder for "D" and "E" cylinders

ANTI-TIPPERS

- Durable steel construction
- Helps prevent users from tipping the wheelchair
- Easy tool-free installation (10260, 10267, & 10268)
- Fits 7/8" and 1" tubing (10271)

URINARY DRAINAGE BAG

- Discreet urological bag holder designed to improve patient dignity
- Constructed with opaque, easy-to-clean fabric
- Durable straps secure to wheelchairs, Geri-chairs, and bedrails