
Ridgewood COLLECTION



COMPOSITE HEAD AND FOOTBOARDS

SPLIT FRAME DESIGN FOR EASY TRANSPORT

CONSTRUCTED FROM DURABLE STEEL



BED FRAME IS ADJUSTABLE FROM 9.8" - 20.5"

EMERGENCY HAND CRANK

EQUIPPED WITH 3" CASTERS THAT ADJUST FROM STANDARD TO LOW HEIGHT

3 MOTOR FULL ELECTRIC HI / LOW HOMECARE BED - 36"

The Dynarex 3-Motor Hi/Low Homecare Bed features a versatile design with pre-mounted motors to allow for optimal convenience during assembly and transport. The bed frame adjusts between standard height and low height to best suit the needs of each individual patient.

- Weight capacity **450 lb.**
- Warranty: Limited lifetime on welds, 5 years on frame, 1 year on all other components

10408	Bed Only	HCPC: E0297	**10408-HR**	W/ Half Rail	HCPC: E0266
10408-FR	W/ Full Rail	HCPC: E0266	**10408-HRFM**	W/ Half Rail / Foam Mattress	HCPC: E0265
10408-FRFM	W/ Full Rail / Foam Mattress	HCPC: E0265	**10408-HRMZ**	W/ Half Rail / Multi-Zone Mattress	HCPC: E0265
10408-FRMZ	W/ Full Rail / Multi-Zone Mattress	HCPC: E0265	**10408-HRTL**	W/ Half Rail / Tri-Layer Mattress	HCPC: E0265
10408-FRTL	W/ Full Rail / Tri-Layer Mattress	HCPC: E0265	**10408-HRIS**	W/ Half Rail / Inner-Spring Mattress	HCPC: E0265
10408-FRIS	W/ Full Rail / Inner-Spring Mattress	HCPC: E0265			

| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 | 1.888 DYNAREX | DYNAREX.COM |