

dynarex®
we care like family™



DYNAREST MULTI-ZONE FOAM PRESSURE MATTRESS

Redistribution Mattress

10420 – 36" × 80"

The DynaRest Multi-Zone Foam Pressure Mattress is constructed from comfortable high-density foam with multiple zones to reduce shear friction, facilitate air flow, and help prevent pressure wounds. Ideal for use at home and in a wide range of healthcare facilities.

- Horizontal precision-cut high-density foam with pressure redistribution zones
- Meets 16 CFR Part 1632 and Part 1633 Standards for Flammability
- Includes removable, fluid-resistant, stretch cover with anti-contamination zipper
- Weight capacity 350 lb.
- 3-year Limited Foam Warranty, 1-year Cover Warranty
- 1 per case
- HCPC: E0184

| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 | 1.888 DYNAREX | DYNAREX.COM |