


dynarex®
we care like family™

CLOSED SYSTEM FOLEY CATHETER TRAYS

ITEMS 5003-5007

Dynarex Closed System Foley Catheter Trays contain the components needed to provide safe, simple catheterization to help minimize tissue irritation and reduce the risk of urinary tract infections.

- Contains essential components for safe and efficient catheterization
- Ready-to-use closed catheter system packed for optimal convenience
- Designed to help reduce the risk of urinary tract infections
- 5004-5007: Includes Foley Catheter (5 cc), 2000 mL Drainage Bag, 10-cc Syringe pre-filled with sterile water, Bed Sheet Clamp, Waterproof Drape, Fenestrated Drape, 2 Medium Powder-Free Vinyl Gloves, 3 PVP Swabsticks in a Pack, Lubricating Jelly, Specimen Container, Patient Label, Packaging Overwrap, Procedure Tray, Insert Label
- 5003: Catheter-free tray contains all components listed above except for the Foley Catheter
- Foley Catheter available in 12 Fr (5004), 14 Fr (5005), 16 Fr (5006), and 18 Fr (5007)
- Sterile, disposable, single-use only for optimal hygiene
- This product contains natural rubber latex
- 10 kits per case


| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 | 1.888 DYNAREX | DYNAREX.COM |



dynarex®
we care like family™

CLOSED SYSTEM FOLEY CATHETER TRAYS

TRAYS INCLUDE



PROCEDURE TRAY



SYRINGE



WATERPROOF DRAPE



FOLEY CATHETER, DRAINAGE BAG & BED SHEET CLAMP



3 PVP SWABSTICKS

FENESTRATED DRAPE



LUBRICATING JELLY



2 VINYL GLOVES



SPECIMEN CONTAINER & PATIENT LABEL



INSERT LABEL

| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 | 1.888 DYNAREX | DYNAREX.COM |