

dynarex®
we care like family™


INSULIN SYRINGES

ITEM 6940-6959

Dynarex Insulin Syringes are high-quality hypodermic syringes designed to effectively deliver insulin. Each disposable syringe features a sterile, stainless-steel needle with a lubricated double-beveled tip for improved patient comfort. The transparent barrel is clearly marked with large-font numbering for accurate dosing and measurement.

- Premium syringe with stainless-steel needle designed to delivery insulin safely and effectively
- Needle is securely attached to help prevent accidental needle pop-off
- Features lubricated thin-walled, double-beveled needle for improved patient comfort
- Transparent barrel with large-font gradations allows for clear visualization
- For use with U-100 insulin
- Non-toxic and non-pyrogenic
- 6940-6949: Individually wrapped, sterile, and disposable for optimal hygiene, 100 syringes per box, 5 boxes per case
- 6951-6959: Convenient multi-pack for quick and efficient dispensing, 10 syringes per bag, 10 bags per box, 5 boxes per case
- Not made with natural rubber latex

| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 | 1.888 DYNAREX | DYNAREX.COM |



dynarex®
we care like family™



INSULIN SYRINGES

U-100

Item	Gauge	Length	Diameter	CC
6940	27 G	½"	0.42 mm	0.5 cc
6945	27 G	½"	0.42 mm	1 cc
6946	28 G	½"	0.37 mm	1 cc
6947	29 G	½"	0.35 mm	1 cc
6948	30 G	½"	0.32 mm	1 cc
6949	31 G	5/16"	0.27 mm	1 cc
6951	28 G	½"	0.37 mm	0.5 cc
6952	29 G	½"	0.35 mm	0.5 cc
6953	30 G	5/16"	0.42 mm	0.5 cc
6954	31 G	5/16"	0.27 mm	0.5 cc
6955	27 G	½"	0.42 mm	1 cc
6956	28 G	½"	0.37 mm	1 cc
6957	29 G	½"	0.35 mm	1 cc
6958	30 G	½"	0.42 mm	1 cc
6959	31 G	5/16"	0.27 mm	1 cc

| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 |
| 1.888 DYNAREX | DYNAREX.COM |