
THE
# MONTCLAIR
COLLECTION

155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 | 1.888 DYNAREX | DYNAREX.COM




## D100 LONG-TERM CARE BED

3 FUNCTION STANDARD BED

The Dynarex D100 Standard Height Bed is constructed from durable steel and features sturdy locking swivel casters, wall bumper guards, and tool-free integrated quick-release head and footboards. A convenient hand pendant allows users to independently adjust bed positioning. This bed is often used in a variety of long-term care facilities including nursing homes, hospitals, rehabilitation centers, and more. Alternate D100 options include additional colors, head and footboards, backup batteries, nurse control panel, and bed rails.

- Weight capacity **450 LB.**
- Warranty: 15 years on all welds, 3 years on electrical and mechanical, and 1 year on all other parts and components

## 12001 BUNDLES

| | AMBER CHERRY BOARDS | MAHOGANY BOARDS | LIGHT OAK BOARDS |
|---|---|---|---|
| **NO RAIL** | 12001-CH | 12001-MH | 12001-OA |
| **COMPOSITE SWING RAIL** | 12001-CHCR | 12001-MHCR | 12001-OACR |
| **METAL SWING RAIL** | 12001-CHMR | 12001-MHMR | 12001-OAMR |
| **PIVOT ASSIST BAR** | 12001-CHPB | 12001-MHPB | 12001-OAPB |
| **ASSIST BAR** | 12001-CHAB | 12001-MHAB | 12001-OAAB |