
THE
MONTCLAIR
COLLECTION

TRUE VERTICAL LIFT FROM 7.8" - 30"

TOOL-FREE REMOVABLE HEAD AND FOOT BOARDS

36" × 80" GRID DECK

EQUIPPED WITH MATTRESS RETAINERS



WALL BUMPER GUARD PROTECTION

INTEGRATED, TOOL-FREE LENGTH EXTENSION 80" - 84" - 88"

EQUIPPED WITH EIGHT 3" TOTAL LOCK SWIVEL CASTERS

D200 LONG-TERM CARE BED

3 FUNCTION LOW BED

The Dynarex D200 Low Bed is a durable, one-piece steel bed frame designed for long-term care use. The frame features true vertical lift that raises up to 30" for increased caregiver convenience and comfort while tending to patients. It can be used with a variety of pressure mattresses to help patients with limited mobility, chronic illness, and pressure wounds. Ideal for use in a wide range of healthcare facilities.

- Weight capacity **450 LB.**
- Quick-release pin lock insert attaches head and footboards with wall bumper guard protection
- Warranty: 15 years on all welds, 3 years on electrical and mechanical, and 1 year on all other parts and components

12002 BUNDLES

	AMBER CHERRY BOARDS	MAHOGANY BOARDS	LIGHT OAK BOARDS	COMPOSITE BOARDS
NO RAIL	12002-CH	12002-MH	12002-OA	N/A
COMPOSITE SWING RAIL	12002-CHCR	12002-MHCR	12002-OACR	12002-CMCR
METAL SWING RAIL	12002-CHMR	12002-MHMR	12002-OAMR	12002-CMMR
PIVOT ASSIST BAR	12002-CHPB	12002-MHPB	12002-OAPB	12002-CMPB
ASSIST BAR	12002-CHAB	12002-MHAB	12002-OAAB	12002-CMAB

155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 | 1.888 DYNAREX | DYNAREX.COM