THE MONTCLAIR COLLECTION

EQUIPPED WITH MATTRESS RETAINERS

COMFORT CHAIR POSITION

INTEGRATED WOODEN HEAD AND FOOT BOARDS

TRUE VERTICAL LIFT FROM 3.9" - 22"

SPLIT FRAME DESIGN FOR TRANSPORT MODE

36" × 80" DECK



AUTO CONTOUR



INCLUDES 10 - FUNCTION HAND CONTROL PENDANT

TRENDELENBURG / REVERSE TRENDELENBURG

EQUIPPED WITH FOUR 3" TOTAL LOCK SWIVEL CASTERS WITH CASTER GUIDE

155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 | 1.888 DYNAREX | DYNAREX.COM

D400 LONG-TERM CARE BED

ULTRA-LOW BED

The Dynarex D400 Ultra-Low Bed is a durable, split-frame bed designed for long-term care use. The bedframe features true vertical lift and a low minimum height to help prevent patient injuries. It can be used with a variety of pressure mattresses to help patients with limited mobility, chronic illness, and pressure wounds. Ideal for use in a wide range of healthcare facilities.

- Weight capacity **450 LB.**
- Constructed from welded steel
- Warranty: 15 years on all welds, 3 years on electrical and mechanical, and 1 year on all other parts and components

12004 BUNDLES

LIGHT OAK BOARDS

NO RAIL	12004-OA
COMPOSITE SWING RAIL	12004-OACR
METAL SWING RAIL	12004-OAMR
PIVOT ASSIST BAR	12004-OAPB
ASSIST BAR	12004-OAAB