

dynarex®
we care like family™

LAB JACKETS WITH POCKETS

ITEMS 2002-2006 & 2012-2016

Dynarex Disposable Lab Jackets are designed to protect skin and clothes from light dirt, splashes, and spills. Constructed from breathable polypropylene material, these durable garments provide a comfortable, non-constricting fit over everyday clothing.

- Offers ANSI/AAMI PB70 Level 1 protection
- Triple-layered polymer material resists moisture, dirt, and spills
- Waist-length garment provides ample upper-body coverage
- Elastic wrists provide a secure fit and reduce the risk of snags and tears
- Buttons in front for reliable closure
- Reuse or dispose as needed
- Not made with natural rubber latex
- Semi-transparent

| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 | 1.888 DYNAREX | DYNAREX.COM |



dynarex®
we care like family™

LAB JACKETS WITH POCKETS

WHITE

2002	Small	10 / bg	3 bgs / cs
2003	Medium	10 / bg	3 bgs / cs
2004	Large	10 / bg	3 bgs / cs
2005	X-Large	10 / bg	3 bgs / cs
2006	2X-Large	10 / bg	3 bgs / cs

BLUE

2012	Small	10 / bg	3 bgs / cs
2013	Medium	10 / bg	3 bgs / cs
2014	Large	10 / bg	3 bgs / cs
2015	X-Large	10 / bg	3 bgs / cs
2016	2X-Large	10 / bg	3 bgs / cs

| 155 CHESTNUT RIDGE ROAD, MONTVALE, NJ 07645 | 1.888 DYNAREX | DYNAREX.COM |