
LARYNGEAL MASK AIRWAYS

Silicone

Dynarex Resp-O_2™ Laryngeal Mask Airways (LMA) are premium-quality respiratory devices designed to help open and maintain airways in unconscious patients. The Laryngeal Mask Airway features a cushioned inflatable cuff with reinforced wire coil to provide a direct, sealed route to the lungs. The Laryngeal Mask Airway is ideal for use by anesthesiologists administering general anesthesia during surgery and other medical procedures. The LMA also provides a simpler alternative to endotracheal intubation in emergency situations, allowing for relatively easy insertion and blind intubation, even by novice operators.

- Inflatable cuff fills oral pharynx to help protect airway
- Soft materials and contoured design help reduce throat irritation
- Ideal for use during the administration of anesthesia
- Allows for relatively easy placement and blind intubation
- Translucent tubing and depth marks allow for accurate placement
- Reinforced with wire coil for optimal support and reduced risk of kinking
- Single-use mask is sterile and x-ray detectable


The **Reinforced** masks feature a **wire coil** for maximum strength, reducing risk of kinking



Inflatable Cuff provides direct sealed routed to lungs

Inflation Pilot Balloon

Valve



1.0 mm Reinforced	1.5 mm Reinforced	2.0 mm Reinforced	2.5 mm Reinforced	3.0 mm Reinforced	4.0 mm Reinforced	5.0 mm Reinforced
REORDER NO. 36285	**REORDER NO. 36286**	**REORDER NO. 36288**	**REORDER NO. 36289**	**REORDER NO. 36290**	**REORDER NO. 36291**	**REORDER NO. 36292**
≤5 KG	**5-10 KG**	**10-20 KG**	**20-30 KG**	**30-50 KG**	**50-70 KG**	**70-100 KG**
STERILE	**STERILE**	**STERILE**	**STERILE**	**STERILE**	**STERILE**	**STERILE**
PACKAGE QUANTITY 5/Case	**PACKAGE QUANTITY** 5/Case	**PACKAGE QUANTITY** 5/Case	**PACKAGE QUANTITY** 5/Case	**PACKAGE QUANTITY** 5/Case	**PACKAGE QUANTITY** 5/Case	**PACKAGE QUANTITY** 5/Case

1.0 mm Non-Reinforced	1.5 mm Non-Reinforced	2.0 mm Non-Reinforced	2.5 mm Non-Reinforced	3.0 mm Non-Reinforced	4.0 mm Non-Reinforced	5.0 mm Non-Reinforced
REORDER NO. 36293	**REORDER NO. 36294**	**REORDER NO. 36295**	**REORDER NO. 36296**	**REORDER NO. 36297**	**REORDER NO. 36298**	**REORDER NO. 36299**
5 KG	**5-10 KG**	**10-20 KG**	**20-30 KG**	**30-50 KG**	**50-70 KG**	**70-100 KG**
STERILE	**STERILE**	**STERILE**	**STERILE**	**STERILE**	**STERILE**	**STERILE**
PACKAGE QUANTITY 10/Case	**PACKAGE QUANTITY** 10/Case	**PACKAGE QUANTITY** 10/Case	**PACKAGE QUANTITY** 10/Case	**PACKAGE QUANTITY** 10/Case	**PACKAGE QUANTITY** 10/Case	**PACKAGE QUANTITY** 10/Case

R250331

dynarex® we care like family™

1.888 DYNAREX | DYNAREX.COM | 155 CHESTNUT RIDGE RD | MONTVALE, NJ 07645

Contact your local sales representative or call to find your nearest distributor.

®2025 Dynarex Corporation