LAB CHOICE

155 Chestnut Ridge Road, Suite 2202
Montvale, NJ 07645
1.888 Dynarex | Dynarex.com





POLYSTYRENE

Item 43050



POLYPROPYLENE

Item 43051

UNIVERSAL CONTAINER

- Designed for safe and efficient collection of both solid and liquid samples
- Constructed from premium polystyrene and polypropylene materials for optimal durability
- Features screw-on cap for a secure fit and reduced risk of leakage
- 30 mL capacity
- Not made with natural rubber latex
- 43050: Maximum centrifugation at 4000×g
- 43051: Maximum centrifugation at 6000×g
- 100 per bag, 4 bags per case

dynarex®
we care like family™